(1) UBIGER

(1) UBIGENE **STR Authentication Report For Cell Line**

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us

Online Inquiry

1. Cell Name: MDA-MB-361

- 2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer. JUBIGEN
- 3. Sample gene genotype analysis result:

	Genotype	e analysis resu	ılts of STR and	Amelogenin l	oci in cells		
STR Loci	STR profile for Sample cell Cell name: STR25032701			STR profile for Cell Bank cell Cell name: MDA-MB-361			
	D5S818	10	11	NBIC	10	11	
D13S317	11	11		11	11	-61	
D7S820	9	12		9	12	- IBIGE	
D16S539	11	12		11	12		
VWA	17	17		17	17		
TH01	9.3	9.3		9.3	9.3		
AMEL	Х	Х		Х	Х		
TPOX	8	11		8	11		
CSF1PO	12	12	UBIC	12	12		
D12S391	21	22	23.0			a E	
FGA	20	24				NIBION	
D2S1338	19	20					
D21S11	30	32.2					
D18S51	12	15					
D8S1179	15	15					

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us Online Inquiry

D3S1358	16	16			
D6S1043	14	18			
PENTAE	8	14	. C	ENE	
D19S433	12.2	14	UBIC		
PENTAD	12	13			-E1
D1S1656	12	15			* I IBIOL

4. Conclusion: The STR results of MDA-MB-361 cells are consistent with the genotypes of MDA-MB-361 cell lines in DSMZ database, the cell ID corresponded to HTB-27, and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com Web: www.ubigene.us



5. Attached Image



Tel: +86 153 6067 3248(Int'l) US Toll free: 855 777 3210 EU Toll free: 800 3272 9252

Korea Toll free: 001 800 3272 9252

Email: info@ubigene.com
Web: www.ubigene.us

Online Inquiry