

STR Authentication Report For Cell Line

1. Cell Name: U373-MG

2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.

3. Sample gene genotype analysis result:

Genotype analysis results of STR and Amelogenin loci in cells						
STR Loci	STR profile for Sample cell			STR profile for Cell Bank cell		
	Cell name: STR24112701			Cell name: U-373 MG		
	Allele1	Allele2	Allele3	Allele1	Allele2	Allele3
D5S818	11	12		11	12	
D13S317	10	11		10	11	
D7S820	10	12		10	12	
D16S539	12	12		12	12	
VWA	16	18		16	18	
TH01	9.3	9.3		9.3	9.3	
AMEL	X	X		X	X	
TPOX	8	8		8	8	
CSF1PO	11	12		11	12	
D12S391	17	22				
FGA	21	25				
D2S1338	22	24				
D21S11	29	29				
D18S51	13	13				
D8S1179	13	15				
D3S1358	16	17				
D6S1043	12	12				

PENTAE	7	10				
D19S433	13	15				
PENTAD	12	12				
D1S1656	16	17.3				

4. Conclusion: The STR results of **U373-MG** cells are consistent with the genotypes of **U373-MG** cell lines in **DSMZ** database, the cell ID corresponded to **HTB-17**, and the STR results **basically matched**. In the test, **Multi allele was found** in this cell line.

5. Attached Image

