

STR Authentication Report For Cell Line

1. Cell Name: SW1573
2. Test Method: DNA was extracted using the genome extraction kit (Axygen), amplified using a 21- STR amplification protocol, the STR loci and gender gene Amelogenin were detected on the ABI 3730XL Genetic Analyzer.
3. Sample gene genotype analysis result:

| Genotype analysis results of STR and Amelogenin loci in cells | | | | | | |
|--|-----------------------------|---------|---------|-------------------------------------|---------|---------|
| STR Loci | STR profile for Sample cell | | | STR profile for Cell Bank cell | | |
| | Cell name: STR25030716 | | | Cell name: SW-1573 [SW-1573,SW1573] | | |
| | Allele1 | Allele2 | Allele3 | Allele1 | Allele2 | Allele3 |
| D5S818 | 12 | 13 | | 12 | 13 | |
| D13S317 | 11 | 11 | | 11 | 11 | |
| D7S820 | 9 | 11 | | 9 | 11 | |
| D16S539 | 12 | 14 | | 12 | 14 | |
| VWA | 20 | 20 | | 20 | 20 | |
| TH01 | 6 | 9.3 | | 6 | 9.3 | |
| AMEL | X | X | | X | X | |
| TPOX | 8 | 11 | | 8 | 11 | |
| CSF1PO | 10 | 12 | | 10 | 12 | |
| D12S391 | 21 | 21 | | | | |
| FGA | 21 | 21 | | | | |
| D2S1338 | 20 | 20 | | | | |
| D21S11 | 29 | 29 | | | | |
| D18S51 | 16 | 16 | | | | |
| D8S1179 | 11 | 11 | | | | |
| D3S1358 | 14 | 18 | | | | |
| D6S1043 | 12 | 17 | | | | |

| | | | | | | |
|---------|----|------|--|--|--|--|
| PENTAE | 5 | 13 | | | | |
| D19S433 | 13 | 13 | | | | |
| PENTAD | 9 | 9 | | | | |
| D1S1656 | 11 | 17.3 | | | | |

4. Conclusion: The STR results of **SW1573** cells are consistent with the genotypes of **SW1573** cell lines in **DSMZ** database, the cell ID corresponded to **CRL-2170**, and the STR results **completely matched**. In the test, **No multi allele was found** in this cell line.

5. Attached Image

